



Frequently Asked Questions: *Canine Parvovirus Antigen Test Kit*

Q: Does vaccination interfere with this test?

A: Vaccination cannot interfere and will not give a false positive result using our test. Professor Ronald Schultz at the University of Wisconsin has performed a study in which he demonstrated that parvovirus shedding in the feces after vaccination is below the level of detection of hemagglutination assay, the gold standard laboratory method for detection of parvovirus.

Q: Why am I having problems getting the test to flow?

A: Most likely, you have used too large of a fecal sample. This makes the mixture in the tube too dense to correctly flow through the test cassette.